

INTRODUCTION TO MDSplus

by
Jeff Schachter

Presented at
General Atomics
San Diego, California

SEPTEMBER 22, 1999



CLASS OUTLINE

Class 1: What is MDSplus and why should you care?

1. MDSplus — the philosophy
2. MDSplus at DIII-D
3. Basic concepts
4. MDSplus expressions
5. What you can do with MDSplus
6. Getting started


Class 2: How to use MDSplus

1. Overview
2. Examining tree structure
3. Reading data out of MDSplus
4. Writing data into MDSplus
5. Designing and building trees

CLASS 2: HOW TO USE MDSplus

- This class presents the nuts and bolts of using MDSplus
- Overview
 - Examining tree (traverser)
 - Opening and closing trees
 - Reading data
 - Using MDSplus in existing tools (ReviewPlus)
 - Writing data to trees
 - Designing and building trees
 - Resources
- Use IDL for MDSplus function calls
 - If you would like C or Fortran information, please e-mail schacht@fusion.gat.com

TRAVERSER MDSplus GUI FOR EXAMINING TREES



USING TCL TO EXAMINE TREES

- **TCL is “Tree Command Language”:** used for browsing and editing trees
- **Like a file system**
 - Current location in tree (“default”)
 - Relative versus absolute paths
 - Set default, dir commands
- **Decompile:** examine contents of node (if not too large!)
- **Show db:** current tree and shot
- **Can do TCL in traverser**
- **Useful for writing macros and generating large trees**
- **Demonstration**

Start TCL from IDL	mdstcl
Open EFIT01 Model	set tree EFIT01/shot=-1/readonly
Directory of tree	dir
Change default	set def .RESULTS.AEQDSK
Examine node contents	decompile WMHD

OPENING AND CLOSING TREES

- Tree must be opened before data within can be accessed
- Opening a tree opens all its subtrees unless you say otherwise
 - Can open subtree directly
- Demonstration

mdsopen, 'EFIT01',97979
- Trees do not need to be closed
 - Most recently opened tree is current context
 - Reopening an already open tree saves overhead
 - Maximum of eight trees open at once, 9th open automatically closes tree on bottom of stack

ASIDE: MDSplus STATUS CODES

- All MDSplus routines in IDL have optional STATUS=status keyword
- If status is odd, success. If status is even, failure
- Use in IDL:

```
mdsopen, tree, shot, /quiet, status=status  
if (status) then begin  
    [ . . . retrieve data . . . ]  
endif else begin  
    print, 'Error opening'+tree+strtrim(shot,2)+mdsgetmsg(status)  
endelse
```
- /quiet keyword suppresses messages
- mdsgetmsg() function returns string explaining error message corresponding to status code

(3) READING DATA

- **Three steps**

- Locate node of interest (tree, path, and shot)
- Open tree
- `mdsvalue(path)`

- **Example: Core Thomson Scattering T_e profile**

Open electrons tree, shot 97979

`mdsopen, 'ELECTRONS', 97979`

Ask for tag `\TSTE_CORE`

`data = mdsvalue('\TSTE_CORE',/quiet,status=status)`

Check status

`print,status+' = '+mdsgetmsg(status)`

NODE REFERENCE

- **Node references can be**

- **Absolute**

- `data = mdsvalue('\ELECTRONS::TOP.TS.BLESSED.CORE:TEMP')`

- **Relative**

- `mdssetdefault,'\ELECTRONS::TOP.TS.BLESSED.CORE'`

- `data = mdsvalue('TEMP')`

- **Tag**

- `data = mdsvalue('\TSTE_CORE')`

- **Recommend that you avoid relative node references for now**

- **Absolute or tag references are easier to understand when examining code**

- **No ambiguity or potential for mistakes**

- `mdssetdefault,'\ELECTRONS::TOP.TS.BLESSED.CORE'`

- `te = mdsvalue('TEMP')`

- `mdssetdefault,'-.DIVERTOR'`

- `t = mdsvalue('DIM_OF(TEMP)')`

- **Gets T_e from core system, but timebase from divertor system**

PASSING PARAMETERS TO MDSVALUE

- `mdsvalue()` function (and `mdsput` procedure) accept parameter arguments
- Primary argument to `mdsvalue()` is a string expression
- Use “\$” within the string to represent a parameter passed in from IDL
- `data = mdsvalue('SOME_FUNCTION($,$)',arg1,arg2)`
- Simple example:
$$\text{arg1} = 2. * !PI * \text{FINDGEN}(101)/100$$
$$\text{data} = \text{mdsvalue}('SIN(\$)',\text{arg1})$$
- Will come up later in use of `GETNCI()` and especially `mdsput`

CACHING AND TDI VARIABLES

- Every time you retrieve a node reference, the tree is read
 - Can be painful if large amount of data
 - O/S does some caching
- Can use TDI variable to “cache” the data in memory
- Example:

No caching	Caching: use TDI variable
te = mdsvalue('\TSTE_CORE')	te = mdsvalue(' _s = \TSTE_CORE ')
t = mdsvalue('DIM_OF(\TSTE_CORE,0)')	t = mdsvalue('DIM_OF(_s,0)')
z = mdsvalue('DIM_OF(\TSTE_CORE,1)')	z = mdsvalue('DIM_OF(_s,1)')
u = mdsvalue('UNITS(\TSTE_CORE)')	u = mdsvalue('UNITS(_s)')
Four separate reads of node	_s is a TDI variable, now equal to the signal \TSTE_CORE

- All TDI variables start with underscore character “_”
 - Distinguish them from node references

USING TDI

- **Math operations**

```
data = mdsvalue(' \TSTE_CORE / 1000. ' )
```

- **Can have expressions containing multiple signals**

- Example: $|I/aB|$ from EFIT

```
mdsopen, 'EFIT01', 97979
```

```
data = mdsvalue(' _s = ABS( \IPMEAS / \AMINOR / \BT0 / 1.e6)')
```

- **Be careful that signals are all the same shape**

- In example, the EFIT signals all have the same timebase

- MDSplus will not prevent you from doing the wrong thing

- **Signal format is lost: no more dimensions**

```
print, mdsvalue('DECOMPILE(_s)')
```

```
Build_With_Units([. . .], "A/m/T")
```

- **But units information retained**

```
print, mdsvalue('UNITS(_s)')
```

```
A/m/T
```

USING GETNCI

- Arguments to GETNCI: GETNCI (node_reference, item_requested, [usage])
- node_reference
 - Some node or set of nodes in a tree
 - Use “\” if absolute path or tag (“\TSTE_CORE”, “\TOP:NAMELIST”)
 - Wildcards:
 - * = all nodes at this level
 - *** = all nodes at this level and below
 - Example: all nodes one level down \TOP.RESULTS.AEQDSK in EFIT tree
“\EFIT01::TOP.RESULTS.AEQDSK:”
 - Example: All RESULTS nodes in EFIT tree
“\EFIT01::TOP.RESULTS***”
 - If using wildcards, return value must be same shape for each node
- item_requested
 - Full name, tag name (minimum path), length, relatives, data
 - Many different possibilities: see TDI help, GETNCI function
- usage
 - If specified, limit wildcard search to nodes of specified type
 - Signal, text, numeric, etc.

USING GETNCI — AN EXAMPLE

- Get names and labels of EFIT signals

- Get node ID numbers (NID) for further reference:

```
nids = mdsvalue('GETNCI("\\EFIT01::TOP.RESULTS.*:*:LABEL","NID_NUMBER"))
```

- Get labels from each node

```
labels = mdsvalue('GETNCI($,"RECORD")',nids)
```

- Get signal names

```
signals = mdsvalue('GETNCI(GETNCI($,"PARENT"),"MINPATH")',nids)
```

➔ Inner GETNCI() gets parent NID numbers

➔ Outer GETNCI() gets minimum path (tag) for each parent NID

GETTING PTDATA FROM MDSplus

- TDI function PTDATA() is interface between MDSplus and PTDATA
- PTDATA() returns signal containing data, timebase, and units
- Example: IP for shot 96021
 - `data = mdsvalue('_s = PTDATA($,$)', "IP", 96021)`
 - `time = mdsvalue('DIM_OF(_s)')`
 - `units = mdsvalue('UNITS(_s)')`
 - ➔ Note parameter passing and caching
- Can also get PTDATA headers
 - `Pthead_IFIX(pointname,shot)`
 - `Pthead_RFIX(pointname,shot)`
 - `Pthead_ASCII(pointname,shot)`
 - `Pthead_INT16(pointname,shot)`
 - `Pthead_INT32(pointname,shot)`
 - `Pthead_REAL32(pointname,shot)`
- Interpreting the headers is up to you!
 - See the DFI web pages <http://lithos.ga.com/dfi/>

(4) USING MDSplus IN EXISTING TOOLS

- **MDSplus is already integrated into GA codes**
 - ReviewPlus
 - EFITtools
 - GAprofiles
 - READA/G
 - IDL “get” routines
 - etc.
- **Plotting pointnames in ReviewPlus and REVIEW**
 - “pointname” — string designating signal of interest
 - MDSplus searched to determine tree and signal
 - If not MDSplus, then pseudopointname or PTDATA
 - ➔ ReviewPlus uses IDL to evaluate data combinations
- **Can also use TDI directly in ReviewPlus**
 - Specify tree name in addition to TDI expression
 - TDI expression evaluator used

(5) WRITING DATA TO TREES

- Practice on “scratch” trees first!
- Basic command: `mdsput, node_reference, expression [. . . arguments]`
- Example:
 - Open scratch tree
`mdsopen,'TEST',10`
 - Put signal into \TOP:SIGNAL
`expression = 'BUILD_SIGNAL(BUILD_WITH_UNITS($,$),*,BUILD_WITH_UNTS($,$))'`
`mdsput,'\TOP:SIGNAL',expression,findgen(100)^2, 'gizmos',findgen(100), 'ms'`
→ Note use of parameters (“\$”) in expression: pass in IDL variables

MORE ON WRITING DATA

- **MDSplus automatically logs username of person writing and data entered**
 - See manually with TCL: `dir/full`
 - Retrieve with GETNCI (“OWNER_ID” and “TIME_INSERTED”)
- **MDSplus does not keep track of previous contents of node**
 - ➔ Be careful not to erase data
- **No built-in mechanism for tracking revisions to data**
 - Thomson scattering stores different revisions
 - But not a generic situation — solved specifically for them
 - Can use elements of their solution for other codes if needed

(6) DESIGNING AND BUILDING TREES

- Get your own data into MDSplus!
- 5 steps:
 - Identify information and data
 - Draw tree structure — how to organize? signal names?
 - Identify dependencies in data
 - Construct tree (using TCL)
 - Write procedure to load it
- Can be iterative process
- Use scratch tree to start
- Once finalized, will become part of D3D tree

CONSIDERATIONS WHEN WRITING DATA

- **Units! Units! Units**
 - Not mandatory, but very helpful
- **Take advantage of indirect referencing for timebases, etc.**
 - Saves space
 - Clearer organization
- **Calibration factors and raw data**
 - Use TDI expressions to allow changing factors without rewriting data
- **Cascading changes**
 - Changing one node affects others
 - Example: “blessed” Thomson revision: \ELECTRONS::TOP.TS:BLESSED_ID
- **Tradeoffs: performance, storage, clarity**
 - In the end — its up to you!

SIGNAL NAMES

- “Pointname” retrieval from MDSplus
 - User does not need to know tree
 - Instead MDSplus function searches all trees for tag
 - Signal tags must therefore be unique across entire DIII-D tree
 - MDSplus only requires unique tags within subtree
- Length limits
 - REVIEW: 10 characters
 - MDSplus node name: 12 characters
 - MDSplus tag: 23 characters
- To be backwards compatible, can only use 10 characters in tag

TCL FOR BUILDING TREES

- **ADD NODE node_reference [/USAGE=usage]**
 - node_reference
 - ★ Can be absolute or relative
 - ★ Use “.NODE” for STRUCTURE nodes
 - usage
 - ★ If node is not STRUCTURE, must specify
 - ★ SIGNAL, NUMERIC, TEXT are most common
- **DELETE NODE node_reference**
 - All nodes under node_reference also deleted
 - Will ask for confirmation
- **ADD TAG node_reference tag_name**
 - tag_name: do not need “\” character
- **WRITE**
 - Save your changes!

TCL TIPS AND TRICKS

- Use macros if adding same group of nodes repeatedly with changes
 - Macro is set of TCL commands grouped together
 - Ask Jeff if interested
- Can write TCL script files
 - Invoke with “@” sign
 - Can pass parameters (VMS syntax)
- Can send TCL commands from IDL
 - Uses mdstcl procedure: `mdstcl,tcl_command_string`
 - Very handy for writing programs to build trees

OTHER TCL CONSIDERATIONS

- **SET NODE command: determine automatic compression and read/write settings**

- General form: SET NODE node_reference /qualifiers

- Qualifiers:

/COMPRESS_ON_PUT : compress data when written

/WRITE_ONCE : cannot rewrite data

/MODEL_WRITE : can write only into MODEL (shot -1) (or /NOMODEL_WRITE)
--

/SHOT_WRITE : can write only into pulse tree (or /NOSHOT_WRITE)

- **Trees can be manually compressed and cleaned**

- Compression: reduce space occupied by data

- Cleaning: remove old data from files

SECURITY AND ACCESS CONTROL

- Right now, all users have all access to all trees
 - BE CAREFUL!
- Will soon move to system where write access is restricted
 - Will be done by group membership
 - Example: all spectroscopists will have write access to SPECTROSCOPY tree
- Access is granted at the subtree level
 - Cannot control access to parts of the tree
- Offsite access will start as read only
 - Can grant write access to specific users from specific sites

(7) RESOURCES

- **DIII-D Data Analysis Group website: <http://fusion.gat.com/comp/analysis>**
 - MDSplus at DIII-D
 - Analysis code documentation
 - Signals documentation
 - Links to MIT site, online MDSplus help
- **E-mail list? Newsgroup?**
 - Wait to see demand
- **Ask Jeff**
 - X4168, 13/413
 - Always happy to answer questions

SUMMARY

- **MDSplus is ready for use!**
- **Little added effort will return many more benefits**
 - Direct access to data
 - Know exactly what you get
 - Retrieve more than just the numbers
- **Get your data into MDSplus**
 - Once in, everyone can see it
 - More efficient for everyone if you do the design work
 - Jeff will help and work on the “big” sets